



BAUSCH+LOMB
STORZ®

Trusted Quality. Leading Design.




CapsuleGuard™ IA




CapsuleGuard™ IA Silicone Irrigation/Aspiration Handpiece

- Flexible Tip
- Two Irrigation Ports
- Flexible Irrigation Tube
- Soft Aspiration Port
- Flexible Sleeve Conforms to Wound to Seal Incision
- One Piece Silicone Construction
- Enhanced Visualization from Semi-transparent Silicone
- No Sharp Edges for Ease of Insertion Into The Eye



CapsuleGuard™ IA



NOW AVAILABLE IN 3 TIP STYLES & 2 INCISION SIZES

Visit www.StorzEye.com/CapsuleGuard for more information and to order online.







EXCELLENT FOR USE IN ALL PHASES OF I/A

"The CapsuleGuard™ has been a game-changer for me in the art of removing cortical material; there is no doubt that, at times, it is more arduous to remove cortex than the nucleus especially with tougher cataracts, smaller pupils, and flimsy posterior capsules. Ever since implementing the CapsuleGuard™ as part of my standard surgical regimen, I feel safer in going after the cortex in all my cases! It's silicone and disposable nature are attributes that I have come to enjoy. Safety and dependability is something we all yearn for and the CapsuleGuard™ does just that."

Jai G. Parekh, MD, MBA, FFAO

Fuqua Alumni Council
The Fuqua School of Business of Duke University,
Class of '08
Managing-Partner, Brar-Parekh Eye Associates
Woodland Park/Edison, NJ

Clinical Associate Professor of Ophthalmology
The New York Eye & Ear Infirmary
New York Medical College, NY
Chief, Cornea & External Diseases
Director, The Research Institute of SJRMC
St. Joseph's Regional Med. Center
Wayne/Paterson, NJ




Stellaris® System

	MICS 1.8mm	Standard Incision 2.2-2.8mm
Straight	85914ST	85912ST
Curved		85913ST
45°	85915ST	85910ST

Other Systems

	MICS 1.8mm	Standard Incision 2.2-2.8mm
Straight	85914S	85912S
Curved		85913S
45°	85915S	85910S



BAUSCH + LOMB
STORZ®

Find the complete range of Bausch + Lomb Storz Ophthalmic Instruments at

www.StorzEye.com

Call Customer Service: 1-800-338-2020