

# Wallace LRI Set

Instrument set provided by R. Bruce Wallace, III, MD of Alexandria, LA

**E0126**

## Wallace LRI Diamond Knife



Single footplate allows easier visibility of knife as it passes through corneal tissue. Handle designed for finger twirling as blade follows arcuate pattern of the limbus. 0.2mm front flat tip. Blade preset at 600 micron. Retractable long handle. Length 120mm, 4.72 inches.

**ET2425**

## Wallace Mendez Degree Gauge



10° interval marks. Designed for surgeons familiar with use of a standard phoropter. Measures 0° to 180° twice in 10° increments. Round handle, length 105 mm, 4.13 inches.

**ET3199**

## Colibri Forceps 0.12mm



0.12mm. 1 x 2 teeth. Serrated handle with polished finish. Overall length: 73mm, 2.9 inches. Manufactured in titanium.

**E7409**

## Plastic Tray - 3 x 6 inches



Size: 3 x 6 inches. Capacity: 4-7 instruments. Depth: 1 inch. Sterilizable plastic tray with silicone matting for supporting and protecting instruments.

